LABTICIAN™
LAB™ AMNIO
 Dehydrated Amniotic Membrane

LabAMNIO™ is a DryTek™ processed, dehydrated, sterilized, human amniotic membrane, prepared in accordance with the Safety of Human Cells, Tissues and Organs for Transplantation Regulations (CTO Regulations), administered by the Biologics and Genetic Therapies Directorate, Health Products and Food Branch, Health Canada.

LabAMNIO™ is a human amniotic membrane allograft that no longer contains living cells. It is processed and sterilized in compliance with appropriate quality management systems to ensure efficacy and safety, and has been declared safe for transplantation. Donor screening, testing, suitability assessment, tissue testing and quarantine are performed under the applicable regulatory requirements as referred in Section 18, 2, 21, and 22 of the CTO Regulations.

Donor selection. Tissue is procured from eligible donors undergoing caesarian section, following voluntary informed consent. All donors are screened using the donor selection criteria as described under the CTO safety requirements in transplantation. The donor's medical history and behavioral risk assessment is obtained and recorded by a qualified medical officer in the hospital. The placenta from the donor is recovered under sterile conditions (operating room), placed into a sterile container and transported under optimal conditions for tissue preservation and protection.

Donor and tissue serological screening. In accordance with the CTO requirements for processing and handling of tissue products, the donor's serological tests include: HUMAN IMMUNODEFICIENCY VIRUS (HIV), HEPATITIS B VIRUS (HBV), HEPATITIS C VIRUS (HCV), and Syphilis/VDRL (Treponemal Specific Assay). Bacteriological testing of the tissue is performed after it is processed and prior to its final sterilisation.

DryTek™ Tissue Processing. LabAMNIO™ is processed in a controlled environment using techniques and protocols designed to prevent contamination and cross-contamination of the tissue. The tissue is carefully washed in a proprietary, specially formulated disinfectant solution to remove any unsolicited biological material. The tissue is dehydrated using a gentle air-dried technique ensuring that its extracellular matrix as well as cellular architecture remain intact. The non-sterile tissue is tested for bioburden prior to sterilization. Once in the desired sizes and thicknesses, the tissue is carefully packaged in an inner and outer pouch, and sterilized to SAL of 10⁻⁶ (Sterility Assurance Level), according to the standards established by the International Organization for Standardization (ISO) and the Association for the Advancement of Medical Instrumentation (AAMI). The double sterile barrier system maintains the sterility of the product and provides protection against biological risks. Each sterilized batch undergoes sterility testing performed by an independent certified laboratory and it released only if the above criteria is met. During the tissue processing, technical quality assurance standards and protocols are rigorously maintained at every step.

Traceability. Recipient records must be maintained for the purpose of tracing tissue post-application in accordance with Sections 43-54 of the CTO Regulations. The patient tracing record is included in the package. The LOT number, Product Code, Donor ID, Tissue size, Manufacturing and Expiry date must be recorded.

Storage. LabAMNIO™ must be stored and maintained in its original package at room temperature until it is ready to be used by the clinician. The product has been validated for storage between 0°C and 40°C.

Shelf-life. Up to 3 years from the date of manufacturing, when stored as recommended.

Ophthalmic usage. LabAMNIO™ is intended for use in ophthalmology as a permanent inlay-graft, temporary onlay-graft or wound covering and, filling-in graft. LabAMNIO™ is an allograft tissue, for homologous use at the recommendation of an eye care professional, for the purpose of supporting the healing processes in various conjunctival and corneal conditions and, to provide a physical barrier in protecting an ocular wound. Depending on the application technique, the tissue can be applied on the ocular surface followed by the placement of a bandage contact lens, or with fibrin glue and/or sutures if necessary.

Recommended instructions for use. These recommendations are designed only to serve as general guidance and they are not intended to replace any of the current institutional protocols, professional clinical judgment and guidelines concerning patient care and management. The use of LabAMNIO™ should be made according to the clinical decision of the physician.

Preparation Instructions

1. Open the box containing the LabAMNIO™.
2. Peel open the outer pouch.
3. Remove the inner pouch using aseptic technique. Both the inner foil pouch and its content may be placed directly onto the sterile field.
4. Peel the inner pouch open only when ready to apply the tissue onto the eye
5. The drier the ocular surface the easier the application.
6. Remove the tissue carefully from the package using a smooth-tip forcep.
7. Place the tissue at the desired ocular area.
8. Gently smoothen out any wrinkles and folds on the tissue using an eye spear.
9. A bandage contact lens is commonly placed over the graft.
Fibrin glue or sutures may be applied if necessary.

Precautions and Contraindications. To avoid any potential complications, LabAMNIO™ should not be used in the presence of an active ocular infection. As with any allograft, complications at the graft site may occur postoperatively. Although donor tissue is evaluated and processed following strict CTO requirements, communicable diseases including those of unknown etiology, transmission of infections agents as well as immune rejection of, or allergic reaction to the tissue are possible.

Adverse events. Adverse reactions or outcomes must be reported immediately to the Labtician Ophthalmics, Quality Assurance and Regulatory Department. Email: qar@labtician.com / Phone: 905-829-0005 ext. 221

Warnings

- Each allograft is intended for single patient use, on a single occasion. DO NOT REUSE.
- Use is limited to eye care professionals only.
- The tissue and packaging should be kept away from the heat.
- Do not freeze the tissue and/or packaging.
- Do not re-sterilize the tissue.
- Do not rinse tissue prior to use.
- Consult instructions for use.
- Discard all open and unused parts of the product.
- Do not use if the package is open, damaged or if mishandling may have caused possible damage and/or contamination.
- Do not use once passed its expiration date. All expired products must be discarded.

The disposal of unused and expired product as well as its packaging must be performed in accordance with accepted medical practice and applicable provincial and federal laws and regulations.

	Caution		Keep Dry
	Do not use if packaging is damaged	 <small>www.Labtician.com/docs</small>	Consult instructions for Use
	Do not reuse		Lot number
	Do not resterilize		Product Code
	Sterilized using gamma		Unique Device Identifier
	Double sterile barrier system		Prescription Only
			Use by date

PRODUCTS	LABTICIAN LABAMNIO 9	LABTICIAN LABAMNIO 9D	LABTICIAN LABAMNIO 12	LABTICIAN LABAMNIO 12V	LABTICIAN LABAMNIO 3x3	LABTICIAN LABAMNIO 3x3D	LABTICIAN LABAMNIO G
	Each LabAMNIO™ product comes packaged in a READYPACK™ kit with a smooth-tip reusable forcep and a box of 10 LabsPEAR™ PVA eye spears						
INDICATIONS	CONJUNCTIVAL SURFACE: Pterygium surgery, chemical burns, conjunctival chalasis, leaking bleb CORNEAL SURFACE RECONSTRUCTION: Non-healing PED neurotrophic ulcers, keratitis sicca, overlay post PRK						
DIMENSIONS	9 mm in diameter 50-micron tissue	9 mm in diameter double thickness 100-micron tissue	12 mm in diameter 50-micron tissue	12 mm in diameter with 2 mm pin-hole 50-micron tissue	3x3cm 50-micron tissue	3x3cm double thickness 100-micron tissue	2x2cm triple thickness 150-micron tissue
ORIENTATION	If the larger projection is distanced clockwise from the smaller one, the epithelial side is facing the cornea.	Double-layered graft. It can be placed in any orientation, as both sides are epithelial surfaces.	If the larger projection is distanced clockwise from the smaller one, the epithelial side is facing the cornea.	If the larger projection is distanced clockwise from the smaller one, the epithelial side is facing the cornea.	If the membrane's cut edge is facing right and away from the surgeon, the stromal side is facing down and the epithelial side is facing up.	Double-layered graft. It can be placed in any orientation, as both sides are stromal surfaces.	Triple-layered graft. Both sides are stromal surfaces.
PRINCIPLES OF SURGERY	<p style="text-align: center;">INLAY OR GRAFT technique OVERLAY OR PATCH technique FILLING-IN OR LAYERED technique</p>						
SHELF LIFE	3 YEARS						
STORAGE CONDITION	Store in a clean, dry environment at ROOM TEMPERATURE						
ADVERSE EFFECTS	The incidence of post-AM transplant microbial infections is very low with DRY AMNIOTIC MEMBRANE						
CONTRA-INDICATIONS	Areas with active or latent infection AM is suitable for transplantation only from an unopened, undamaged package AM is intended for SINGLE-USE ONLY DO NOT RE-STERILIZE						

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Usage ophtalmique LabAMNIO^{MC} is intended for use in ophthalmology as a permanent inlay-graft, temporary onlay-graft or wound covering and, filling-in graft. LabAMNIO^{MC} is an allograft tissue, for homologous use at the recommendation of an eye care professional, for the purpose of supporting the healing processes in various conjunctival and corneal conditions and, to provide a physical barrier in protecting an ocular wound. Depending on the application technique, the tissue can be applied on the ocular surface followed by the placement of a bandage contact lens, or with fibrin glue and/or sutures if necessary.

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Avançons la vision

Labtician Ophthalmics, Inc.,
2140, Winston Park Drive, Suite 6,
Oakville (Ontario) L6H 5V5 CANADA
Health Canada CTO #: 100040

Source Establishment:
Biogentek Lifescience Pvt Ltd
517 Phase 9 Industrial Area, Sector 66 Mohali, Punjab, India, 160062
Health Canada CTO #: 100268