The AHMED® ClearPath Glaucoma Drainage Device

Model CP250 and CP350

NEW WORLD
MEDICAL

CE
0459

Rx_{only} **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

Index

English (en)	4
Español (es)	8
Français (fr)	12
Italiano (it)	16
Nederlands (nl)	20
Português (pt)	24

	ENGLISH	ESPAÑOL	FRANÇAIS	ITALIANO	NEDERLANDS	PORTUGUÊS
	Manufacturer	Fabricante	Fabricant	Produttore	Fabrikant	Fabricante
	Date of manufacture	Fecha de fabricación	Date de fabrication	Data di produzione	Fabricagedatum	Data de fabrico
	Use-by date	Fecha de caducidad	Utiliser avant le	Data di scadenza	Uiterste gebruiksdatum	Prazo de validade
	Serial number	Número de serie	Numéro de série	Numero di serie	Serienummer	Número de série
	Batch Code	Código de lote	Code du lot	Codice lotto	Partijnummer	Código do lote
	Authorized representative in the European Community	Representante autorizado en la Comunidad Europea	Représentant autorisé dans la Communauté européenne	Rappresentante autorizzato nella Comunità Europea	Gemachtigd vertegenwoordiger in de Europese Gemeenschap	Representante autorizado na Comunidade Europeia
	Consult instructions for use	Consultar las instrucciones de uso	Consulter le mode d'emploi	Consultare le istruzioni per l'uso	Raadpleeg de gebruiksaanwijzing	Consultar as Instruções de utilização
	Read all warnings and precautions in instructions for use	Leer todas las advertencias y precauciones incluidas en las instrucciones de uso	Lire tous les avertissements et toutes les mises en garde figurant dans le mode d'emploi	Leggere tutte le avvertenze e precauzioni nelle istruzioni per l'uso	Lees alle waarschuwingen en voorzorgsmaatregelen in de gebruiksaanwijzing	Leia todos os avisos e precauções nas instruções de utilização
	Do not reuse	No reutilizar	Ne pas réutiliser	Non riutilizzare	Niet hergebruiken	Não reutilizar
	Do not re-sterilize	No reesterilizar	Ne pas restériliser	Non risterilizzare	Niet opnieuw steriliseren	Não reesterilizar
	Sterilized using irradiation	Esterilizado mediante radiación	Stérilisé par rayonnement	Sterilizzato per irraggiamento	Gesteriliseerd met behulp van straling	Esterilizado com irradiação
	Do not use if package is damaged	No utilizar si el envase está dañado	Ne pas utiliser si l'emballage est endommagé	Non utilizzare se la confezione è danneggiata	Niet gebruiken als de verpakking beschadigd is	Não utilizar se a embalagem estiver danificada
	MR safe	MR Safe (seguro con la resonancia magnética)	Compatible avec l'IRM	Compatibile con la risonanza magnetica	MRI-veilig	Segurança na presença de RMN
	Upper limit of temperature	Límite superior de temperatura	Límite supérieure de température	Límite superiore di temperatura	Bovengrens voor temperatuur	Límite superior de temperatura
	Catalogue Number	Número de catálogo	Référence catalogue	Numero di catalogo	Catalogusnummer	Número de catálogo

Figure 1 / Figura 1 / Figuur 1


Figure 2 / Figura 2 / Figuur 2

3.

Figure 3 / Figura 3 / Figuur 3


Figure 4 / Figura 4 / Figuur 4


Figure 5 / Figura 5 / Figuur 5



Figure 6 / Figura 6 / Figuur 6

4.


The AHMED® ClearPath Glaucoma Drainage Device

Device Description:


The AHMED® ClearPath glaucoma drainage device is a non-valved drainage device designed to shunt aqueous in eyes suffering from refractory glaucoma. The implant consists of a medical grade silicone tube secured near the anterior suture points of the medical grade silicone episcleral plate. Two models are available that cover surface areas of approximately 250 mm² and 350 mm². The anterior suture points are located on the anterior side of the plate, flanking each side of the tubing track. The silicone plate is barium impregnated to increase ultrasound resolution and identification with CT scans, MRIs and plain skull films. The plate conforms to the shape of the globe at its equator and provides a surface from which fluid can be dispersed. The device is for single use only.

Each AHMED® ClearPath glaucoma drainage device is supplied with a 23-gauge hypodermic needle and a 4-inch 4-0 polypropylene ripcord, for optional use, in a sterile, sealed double-pouch. The use or non-use of both the 23-gauge hypodermic needle and/or the ripcord do not alter the standard care for the implantation of the AHMED® ClearPath glaucoma drainage device.

Model CP250

Model CP350

Model	Surface Area*	Plate Length	Plate Width (Max)	Plate Thickness (Max)	Tube Length	Tube ID	Tube OD
CP250	250 mm ²	16.70 mm	18.47 mm	0.86 mm	32 mm	0.305 mm	0.635 mm
CP350	350 mm ²	16.48 mm	30.47 mm	0.86 mm	32 mm	0.305 mm	0.635 mm

* Surface area does not include the area covered by the suture arms.

Indications for Use:

The AHMED® ClearPath glaucoma drainage devices are indicated for the management of refractory glaucoma where previous surgical treatment has failed or is not expected to provide satisfactory results. Such refractory glaucoma may include but is not limited to: neovascular glaucoma, congenital or infantile glaucoma, and refractory glaucoma resulting from aphakia or uveitis.

Contraindications:

Bacterial conjunctivitis, bacterial corneal ulcers, endophthalmitis, no light perception, orbital cellulitis, bacteremia or septicemia and active scleritis.

Possible Complications and Adverse Events:

Serious or hemorrhagic choroidal detachments, hyphema, hypotony, flat or shallow anterior chamber, phthisis bulbi, retinal detachment, endophthalmitis, conjunctival erosion, wound leak, tube contacting cornea, tube block by iris or vitreous, tube migration anteriorly or posteriorly, corneal edema, need for secondary surgical intervention, uveitis and diplopia.

Warnings & Precautions:

1. Do not use the device if the sterile package integrity has been compromised.
2. Do not re-sterilize the implant.
3. Tampering with the implant can cause malfunctions.
4. Do not use the 4-inch 4-0 polypropylene ripcord as a suture.
5. The AHMED® ClearPath glaucoma drainage devices are intended for single use only.
6. Use of surgical tools that can/may cause damage to the device are not recommended.
7. Do not store at temperatures above 45°C (113°F).

How Supplied/Expiration Date:

Each AHMED® ClearPath glaucoma drainage device is supplied with a 23-gauge hypodermic needle and a 4-inch 4-0 polypropylene ripcord (preloaded in the tube lumen), for optional use, in a sterile, sealed double-pouch. Both pouches are held within a product specific cardboard box. Product Information, Patient Chart Labels, Patient ID card and a New World Medical Implant Notification Card are also enclosed within the sterile package. The device has been terminally sterilized by gamma radiation. Sterility is assured provided that the peel pouch has not been opened or damaged and the sterility expiration date has not lapsed. The peel pouch is intended to be opened using a sterile technique, allowing the implant to be placed into the sterile field. Expiration date is indicated on the outside of the box. The implant

should not be used after the expiration date indicated. **Note: The manufacturer disclaims all warranties expressed or implied, including but not limited to suitability for a particular purpose.**

Directions for Use:

Preparation:

1. The surgeon should examine the implant prior to use for integrity, proper size, model, and expiration date.
2. Open the sterile product pouch and remove the device in a sterile environment using a sterile technique.
3. Examine the implant to ensure that particles are not present, as it can become electrostatically charged upon opening the package. Rinse implant in sterile saline, if required.

Implantation:

The surgeon should have adequate training and familiarity with the implantation of glaucoma drainage devices, generally, and the Ahmed family of devices, specifically, prior to implantation of the AHMED® ClearPath glaucoma drainage device.

1. (**Fig. 1 Optional Ripcord**) The tube must be ligated near the tube-plate junction preferably with a 7-0 vicryl dissolvable suture to occlude the tube. One technique to occlude the tube involves the placement of a releasable ripcord in the tube lumen at the time of the initial surgery. A 7-0 vicryl dissolvable suture is placed around the tube (and consequently around the ripcord) and tightened to prevent any flow through the tube and around the ripcord. A 4-inch 4-0 polypropylene ripcord (preloaded in the tube lumen) is provided with the implant.
2. (**Fig. 2**) A conjunctival peritomy incision is made. A pocket is formed superiorly, avoiding the superior rectus and oblique muscles, with blunt dissection of Tenon's capsule from the episclera.
3. (**Fig. 3**) The implant is placed away from the limbus under the conjunctiva and tenon's capsule between the rectus muscles (the endplate of model CP350 is wider than the model CP250 endplate and has portions that fit under adjacent rectus muscles) and sutured to the episclera. The leading edge of the endplate should be ~ 8-10 mm from the limbus resulting in the suture arms being 6-8 mm from the limbus.

6.

Directions for Use: (continued)

- (Fig. 4)** The drainage tube is trimmed to permit a 2-3 mm insertion of the tube into the anterior chamber. The tube should be bevel cut to an anterior angle of 30° to facilitate insertion.
- (Fig. 5)** The AC is entered at 1-2 mm way from the limbus with a sharp 23-gauge needle to create a needle track, parallel to the iris. In addition, a paracentesis should be performed to allow for quick AC reformation in case of a flattened anterior chamber.

Caution: Care must be taken to ensure that the drainage tube does not contact the iris or corneal endothelium after insertion.
- (Fig. 6)** The drainage tube is inserted approximately 2-3 mm into the anterior chamber through the previously created needle track.
- After surgery, the body's wound healing mechanism results in the formation of a fibrous capsule (bleb) composed of scar tissue around the implant, this scar tissue resists fluid outflow. In approximately 4-6 weeks the 7-0 vicryl suture dissolves. At that time, under a magnified view of a slit lamp biomicroscope or operating microscope the conjunctiva can be entered over the temporal side of the ripcord and the ripcord be removed using jeweler's forceps. Following the ripcord removal, the IOP should be checked and slit lamp and fundoscopic examination should be performed.

Explant:

If removal of the device is necessary, standard surgical technique can be used to isolate the plate and tube from surrounding conjunctiva and Tenon's capsule.

- Any adhesions between the device and surrounding tissues should be amputated.
- Pre-existing sutures are cut and removed to allow for free mobility of the plate and tube.
- The silicone tube is then grasped and removed from the sclerotomy tunnel.
- The entire device is then removed from the operative site.
- The sclerotomy tunnel is sutured closed and water tight followed by re-approximation of the conjunctiva and Tenon's capsule using the surgeon's preferred suture and technique.

Return of Damaged/Unsterile Product:

Contact New World Medical (NWM) for product return policy. Return the device with proper identification and reason for return. Label the return as a biohazard.

Contact Number: (909) 466-4304

Return Address: New World Medical
10763 Edison Court
Rancho Cucamonga
California 91730, USA

Patient Information:

It is recommended that each patient receive proper information regarding glaucoma implants and/or glaucoma drainage devices prior to the decision to implant is made. The implanting surgeon should describe to the patient all the risks involved with glaucoma drainage device surgery. At the end of surgery, the Implant Card should be provided to the patient.

Reporting:

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as product related and that were not previously expected in nature, severity or incidence must be reported to NWM. Report any U.S. product related event to the toll free NWM number at (800) 832-5327. For countries outside of the U.S. please contact your local company representative.

Disclaimer of Liability:

NWM is not responsible for any technique(s) performed which are not stated in the instructions for use to implant the product.

Clinical Data:

Although the AHMED® ClearPath has not been clinically evaluated for performance, the following reference in the published clinical literature summarizes the performance of the predicate device Baerveldt Glaucoma Implant.

- o *Budenz, Donald L., et al. "Five-year treatment outcomes in the Ahmed Baerveldt comparison study." Ophthalmology 122.2 (2015): 308-316.*

The following additional clinical literature is referenced in support of differences in design of the AHMED® ClearPath as compared to the predicate Baerveldt Glaucoma Implant.

- Placement of Device: This study was conducted to determine optimal positioning of various glaucoma drainage devices posterior to the limbus while maintaining a safe distance relative to the optic nerve.
 - o *Kahook, Malik Y., et al. "Location of glaucoma drainage devices relative to the optic nerve." British Journal of Ophthalmology 90.8 (2006): 1010-1013.*
- Implant Design – Ridge: The following studies evaluated the safety and effectiveness of the AGV-FP7 reference device to the predicate Baerveldt Glaucoma Implant and indicate that a lack of a ridge does not significantly impact implant performance.
 - o *Budenz, Donald L., et al. "Five-year treatment outcomes in the Ahmed Baerveldt comparison study." Ophthalmology 122.2 (2015): 308-316.*
 - o *Budenz, Donald L., et al. "Postoperative complications in the Ahmed Baerveldt comparison study during five years of follow-up." American Journal of Ophthalmology 163 (2016): 75-82.*
- Use of Ripcord for Tube Occlusion: A ripcord for optional use is provided within the tube lumen of the Ahmed ClearPath device. Tube occlusion via intraluminal ripcord is a common technique to prevent early postoperative hypotony in non-valved glaucoma drainage devices. Below are examples of medical publications that discuss this technique, its effectiveness and its impact on implant performance.
 - o *Sherwood, Mark B., and M. Fran Smith. "Prevention of early hypotony associated with Molteno implants by a new occluding stent technique." Ophthalmology 100.1 (1993): 85-90.*
 - o *An, Selena J., et al. "Scheduled postoperative ripcord removal in Baerveldt 350 implants corresponds with increased complications: a prospective, randomized trial." Investigative Ophthalmology & Visual Science 59.9 (2018): 2073-2073.*

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o *Kahook, Malik Y., et al. "Location of glaucoma drainage devices relative to the optic nerve." British Journal of Ophthalmology 90.8 (2006): 1010-1013.*

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o *Budenz, Donald L., et al. "Five-year treatment outcomes in the Ahmed Baerveldt comparison study." Ophthalmology 122.2 (2015): 308-316.*

o *Budenz, Donald L., et al. "Postoperative complications in the Ahmed Baerveldt comparison study during five years of follow-up." American Journal of Ophthalmology 163 (2016): 75-82.*

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o *Sherwood, Mark B., and M. Fran Smith. "Prevention of early hypotony associated with Molteno implants by a new occluding stent technique." Ophthalmology 100.1 (1993): 85-90.*

o *An, Selena J., et al. "Scheduled postoperative ripcord removal in Baerveldt 350 implants corresponds with increased complications: a prospective, randomized trial." Investigative Ophthalmology & Visual Science 59.9 (2018): 2073-2073.*

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