F4H5[®] WashOut

Deutsch	1
English	2
Français	3
Italiano	4
Español	5
Português	6
Ελληνικά	7
Nederlands	8
Dansk	9
Norsk	10
Svenska	11
Suomi	12
Eesti keel	13
Latviešu	14
Lietuviškai	15
Islenska	16
Polski	17
Česky	18
Slovenčina	19
Magyar	20
Slovensko	21
Română	22
Български	23
Русский	24

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Vial

Spritze

F4H5® WashOut is a biocompatible solvent for removing silicone oil residue from the retina and for cleaning intraocular lenses after silicone oil tamponades.

Application and posology

After removal of a silicone oil tamponade F4H5® is injected directly into the eye. Standard needles can be used for injection. The quantity needed is determined by the surgeon. The F4H5® WashOut solution must be brought into contact with the silicone oil drops to be dissolved in order to ensure removal. At the end of surgery F4H5® must be completely removed from the eye and, if necessary, replaced by a fluid for long-term tamponade.

Composition

100 % perfluorobutylpentane [C₄F₉C₅H₁₁], ultrapure

Physical properties

Molecular weight	[g/mol]	290.17
Boiling point	[°C]	152
Density	[g/cm ³ at 25 °C]	1.28
Refractive index	[at 20°C]	1.320
Surface tension	[mN/m at 25 °C]	18.3
Interfacial tension	[mN/m at 25 °C]	43.03
Viscosity	[mPas at 25 °C]	1.05

Other Remarks

F4H5® remaining in the vitreous cavity can be completely removed using a vacuum instrument.

Contraindications

F4H5® WashOut must not be used in eyes with severe injuries because the product may leak into the orbital cavity. F4H5® WashOut must not be used in eyes with uveitis because these eyes tend to develop more severe inflammatory conditions. Side effects are unknown.

Interactions

F4H5® is completely soluble in silicone oils. This allows effective removal of large quantities of silicone oil. Interactions with other agents are not known.

Side effects

F4H5® is well tolerated inside the eye. No side effects are known to date.

Form of administration

F4H5® WashOut is available in 5 mL units.

Storage conditions

Store F4H5® WashOut at room temperature.

Sterilisation

Heat sterilisation with subsequent autoclaving.

Cautions

Do not use this product after the expiry date. The expiry date is shown on the vials and on the packaging. Only use undamaged packages. Vials are only intended for single use on one patient. Repeated use can lead to microbial contamination, which may cause a severe infection in the patient. F4H5® WashOut should be injected with the sterile disposable syringe enclosed. If other sterile disposable syringes are used, make sure F4H5® is injected immediately after filling the syringe. Otherwise the syringe may become awkward to use, making smooth delivery difficult. The reason for this is that F4H5® dissolves the silicone coating on conventional disposable syringes.

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FLUORON GmbH