

# Accu-Temp<sup>®</sup> Cautery

Specialty Cauterization

Accu-Temp<sup>®</sup> Cauteries provide a broad selection of tip styles for cauterization in a wide variety of ophthalmic procedures.



For ordering information call:  
1.866.906.8080/1.781.906.8080  
or visit [www.beaver-visitec.com](http://www.beaver-visitec.com)

 **Beaver Visitec**  
*Keeping Your Vision in Sight*

# Accu-Temp® Cautery

## Specialty Cauterization

### High Temperature Cautery



Description	Item #
1205°C (2200°F), 1/2" shaft with fine tip, 10/bx	8442000
1205°C (2200°F), 1/2" shaft with loop tip, 10/bx	8443000
1205°C (2200°F), 2" shaft with loop tip, 10/bx	8444000
1205°C (2200°F), 2" shaft with fine tip, 10/bx	8445000
982°C (1800°F), 1/2" shaft with elongated Vasector fine tip, 10/bx	8441500

### Low Temperature Cautery



Description	Item #
677°C (1200°F), 1/2" shaft with fine tip, 10/Bx	8441000

### Variable Temperature Cautery - Dial Controlled



Description	Item #
454°-649°C (849°-1200°F), 1/2" shaft with fine tip, 10/bx	8440000
813°-1148°C (1495°-2098°F), 1/2" shaft with fine tip, 10/bx	8446000

### Interchangeable Tip Cautery - Reusable Handles, Single-Use Tips



Description	Item #
Variable low temperature handle (requires single-use tips), 10/Bx	8461000
Fixed high temperature handle (requires single-use tips), 10/Bx	8461020
Single-use fine tip, 25/Bx	8461050
Single-use loop tip, 25/Bx	8461060



**Beaver-Visitec International, Inc.**  
 411 Waverley Oaks Road, Waltham, MA 02452  
 Tel: 1.866.906.8080 Fax: 1.866.906.4304  
[www.beaver-visitec.com](http://www.beaver-visitec.com)