










Visitec® Single-Use Microsurgical Instruments

Meeting the Needs of Ophthalmic Surgery

Description	Item #	Description	Item #	Description	Item #
 <p>Vannas Scissors Curved, 9.5 mm blades with sharp pointed tips. 90 mm long. 10 per box.</p>	581422US	 <p>Corneal Scissors [Castroviejo] LEFT. 9 mm blades, blunt tips. 100 mm long. 10 per box.</p>	581469US	 <p>Westcott Scissors 23 mm blades with blunt tips. 105 mm long. 10 per box.</p>	581437US
 <p>Corneal Scissors [Castroviejo] RIGHT. 9 mm blades, blunt tips. 100 mm long. 10 per box.</p>	581470US	 <p>Westcott Scissors 16 mm blades with blunt tips. 112 mm long. 10 per box.</p>	581472US	 <p>Strabismus Scissors Curved, 25 mm blades. 104 mm long. 10 per box.</p>	581471US
 <p>Chalazion Forceps [Desmarres] Oval solid 26 mm wide lower plate. Upper plate oval with inside dimensions 15 mm x 23 mm. 105 mm long. 10 per box.</p>	581450US	 <p>Chalazion Curette Double ended, sharp, 1.50 mm and 2.50 mm diameter. 130 mm long. 10 per box.</p>	581461US	 <p>Strabismus Hook [Graefe] Straight blunt tip. 135 mm long. 10 per box.</p>	581479US



Customized Single-Use Instrument Tray

Many Beaver® and Visitec® single-use microsurgical knives, instruments and cannulae are available in a configurable, customized instrument tray. Customized Single-Use Instrument Trays are designed to:

- Increase efficiency: instruments, knives and cannulae are packaged in a sterile, ready to use tray.
- Eliminate the need for costly and time consuming re-sterilization.
- Simplify ordering and reduce inventory: stocking of individual products is eliminated.



Beaver-Visitec International, Inc.

411 Waverley Oaks Road, Waltham, MA 02452

Tel: 1.866.906.8080 Fax: 1.866.906.4304

www.beaver-visitec.com

BVI, BVI Logo and all other trademarks (unless noted otherwise) are property of a Beaver-Visitec International ("BVI") company © 2011 BVI

590024-0311

Visitec® Single-Use Microsurgical Instruments

Single Source Convenience

Beaver-Visitec is proud to offer a wide variety of stainless steel instruments combining the precision you expect with the convenience of single-use.













For ordering information call:
1.866.906.8080/1.781.906.8080
or visit www.beaver-visitec.com

bvi Beaver Visitec
Keeping Your Vision in Sight

Visitec® Single-Use Microsurgical Instruments

Visitec® single-use instruments are designed to reduce the risk of infection transmission and help you comply with today's safety requirements and recommendations.

Description	Item #	Description	Item #	Description	Item #
 <p>Eye Speculum Open wire loops. 40 x 30 mm, 14 mm blades. 0.8 mm wire size. 10 per box.</p>	581413	 <p>Eye Speculum [Barraquer] Adult, solid 13 mm blades. Width 31 mm. 10 per box.</p>	581474US	 <p>Capsulorhexis Forceps [Utrata] Straight, cystotome tip, 13 mm long bend, flat handle. 95 mm long. 10 per box.</p>	581451US
 <p>Eye Speculum Open wire loops. 40 x 30 mm, 14 mm blades. 0.9 mm wire size. 10 per box.</p>	581443US	 <p>Eye Speculum [Lancaster] Solid blades with spring controlled spread. Blades 15 mm. 80 mm long. 10 per box.</p>	581473US	 <p>Capsulorhexis Forceps [Masket] Curved, cystotome tip, bend 11 mm long, flat handle. Designed to fit in 1.8 mm incision. 100 mm long. 10 per box.</p>	581482US
 <p>Eye Speculum [Barraquer] Adult, closed loop wire 14 mm blades. Width 28 mm. 10 per box.</p>	581477US	 <p>Eye Speculum [Castroviejo] Fenestrated blades with screw controlled spread. 84 mm long. 10 per box.</p>	581475US	 <p>Capsulorhexis Forceps [Masket] Curved, cystotome tip, 11 mm long, round handle. Designed to fit in 1.8 mm incision. 108 mm long. 10 per box.</p>	581484US
 <p>Eye Speculum [Kratz/Barraquer] Open wire loops, 40 x 26.5 mm. Retraction size 13 x 4.5 mm. 3 per box.</p>	585094	 <p>Eye Speculum [Lieberman] Open 14 mm wire blades with screw controlled spread. 10 per box.</p>	581478US	 <p>Capsulorhexis Forceps [Utrata] Straight, cystotome tip, 11 mm long, flat handle. Designed to fit in 1.8 mm incision. 100 mm long. 10 per box.</p>	581481US
 <p>Eye Speculum [Barraquer] Pediatric, closed loop wire 10 mm blades. Width 25 mm. 10 per box.</p>	581476US	 <p>Bonn Forceps 0.12 mm 1 x 2 teeth with 5 mm platform, flat handle. 100 mm long. 10 per box.</p>	581409US	 <p>Capsulorhexis Forceps [Utrata] Straight, cystotome tip, 11 mm long, round handle. Designed to fit in 1.8 mm incision. 108 mm long. 10 per box.</p>	581483US
 <p>Eye Speculum [Wortham] Pediatric, closed loop, 34 x 21 mm. Retraction size 10 x 3.5 mm. 3 per box.</p>	585203	 <p>Capsulorhexis Forceps Curved, cystotome tip, 13 mm long bend, flat handle. 100 mm long. 10 per box.</p>	581410US	 <p>Colibri Forceps [Barraquer] 0.12 mm 1 x 2 teeth. 77 mm long. 10 per box.</p>	581456US

Designed with safety in mind:

- A brand new, sterile instrument for each case ensures patient safety
- Precision stainless steel instruments with the convenience of single-use
- Easy to use, 100% traceability, eliminates the need for costly resterilization

Description	Item #	Description	Item #	Description	Item #
 <p>Colibri Forceps [Castroviejo] 0.12 mm 1 x 2 very delicate teeth, flat handle. 108 mm long. 10 per box.</p>	581457US	 <p>Moorfield Forceps Straight with 14 mm grooved section, flat handle. 120 mm long. 10 per box.</p>	581436US	 <p>Nucleus Manipulator [Drysdale] Flattened paddle tip, 1.10 mm tip. 115 mm long. 10 per box.</p>	581465US
 <p>Jewelers Forceps No. 3. Straight forceps. 100 mm long. 10 per box.</p>	581458US	 <p>Vitreoretinal Forceps Straight, platform on tips. Tip length 31 mm. Designed to fit 23G incision. 146 mm long. 5 per box.</p>	581480US	 <p>Nucleus Rotator [Pisacano] Angled, 0.70 mm diameter flat, blunt tip. 115 mm long. 10 per box.</p>	581466US
 <p>Tying Forceps Curved, 6.0 mm tying platform, flat handle. 96 mm long. 10 per box.</p>	581411US	 <p>Phaco Chopper [Nagahara] 1.4 mm tip with blade offset at 45 degree angle. 115 mm long. 10 per box.</p>	581467US	 <p>Capsule Retractor [Mann/Hunkeler] Angled, 0.6 mm diameter semi-ball. 115 mm long. 10 per box.</p>	581462US
 <p>Micro Tying Forceps Straight, 6.0 mm tying platform, flat handle. 105 mm long. 10 per box.</p>	581459US	 <p>Burrato's Hook Angled, 0.25 mm diameter, blunt tip. 120 mm long. 10 per box.</p>	581412US	 <p>Needle Holder Curved 6 mm jaws. 105 mm long. WITHOUT LOCK. 10 per box.</p>	581417US
 <p>Suture Forceps [Paufigue] Straight shaft with 1 x 2 teeth, 6.0 mm tying platform, flat handle. 85 mm long. 10 per box.</p>	581460US	 <p>Straight Micro Manipulator Angled, 0.35 mm diameter, blunt tip. 120 mm long. 10 per box.</p>	581423US	 <p>Needle Holder [Barraquer] Curved 6 mm jaws. 120 mm long. WITH LOCK. 10 per box.</p>	581468US
 <p>MacPherson Forceps Angled 12 mm from tip to bend, flat handle. 100 mm long. 10 per box.</p>	581416US	 <p>Lens Manipulator [Lester] Angled, 0.44 mm diameter hourglass shape tip. 115 mm long. 10 per box.</p>	581464US		