

The Macula is the central 10% of the Retina, and is the area responsible for the sharpest vision. Through natural degeneration, a hole can form in the Macula causing a decrease in vision and blank spots or blurry spots. In order to try to repair the damage, an Ophthalmologist may have to perform Vitrectomy surgery.

During surgery, the clear vitreous is removed from the center of the eye.

The eye is then filled with a long acting gas bubble that "floats" in the eye and (when positioned correctly) presses against the hole, flattening it on the retina and keeping it dry. The gas bubble slowly dissolves and is replaced naturally by ocular fluid.

Since the gas bubble floats and moves freely to the highest point in the eye, it is very important that you follow the instructions you have been given regarding head positioning following surgery. Keeping the macula dry and positioned properly during the healing process is extremely important.

FACE DOWN

POSITIONING DEVICES

How you can make your retinal surgery a success

LABTICIAN
ophthalmics

Foam Positioning System PRODUCT CODE - AZ701005

Our Foam Positioning System provides ergonomic support for the body when used on the bed or floor by allowing for a more normal curve to the lower back. Proper face down positioning promotes healing while easing stress on your neck, shoulders, spine and lower back.

- Head and face cradle provides generous breathing channel and can be used separately on a table.

The Mobile Support PRODUCT CODE - K501

The Mobile Support System allows you limited mobility and reminds you to stay face down. This product is intended for temporary support only. It will reduce neck strain and assist in maintaining the recommended face down position by resting the weight of your head on a foam support.

- Reduces neck strain
- Adjustable
- Affordable

Head Positioning (Prism) Glasses

PRODUCT CODE - 8080

These glasses allow you to properly position your head while allowing you to watch television, see faces, etc.

- Lightweight, precision ground prisms provide a clear, sharp image.
- Tortoiseshell type plastic frames.
- Removable head strap with hook and loop closure to hold glasses in place.
- Can be used over most regular glasses.

Combo Vitrectomy Kit PRODUCT CODE - KCVK-2

COMBINATION PRODUCT PURCHASE
(Our most popular choice)

Foam Positioning System, Mobile Support Device and Prism Glasses

To order call Toll Free
1-800-265-8391

Labtician Ophthalmics, Inc.
2140 Winston Park Drive,
Unit 6 Oakville, ON,
L6H 5V5 Canada
Tel: 905.829.0055
Fax: 905.829.0056

info@labtician.com
www.labtician.com

VISA & MASTERCARD

LABTICIAN
ophthalmics

PHYSICIAN
PARTICIPATION
CODE OF
CONDUCT

