


Their Advantages Are Easy To Grasp

Broader contact with pupil margin




The MST Iris Hooks Catalog number is:

- MIH-0001
- They are provided with five hooks per package
- Six packages per box

MST Iris Hooks are designed to enter through a 0.6mm paracentesis.

The MST Iris Hooks are best controlled by holding them on top and bottom as opposed to their sides. Holding the top and bottom is recommended to avoid twisting.

Cross Section of Strands



Hold top and bottom not sides

 *Iris Hooks*

For more information: 1.888.279.3323 info@microsurgical.com

Patent Pending