



Other Products

last updated 11.13



DDMS™ Diamond Dusted Membrane Scrapers

The Original Tano Diamond Dusted Membrane Scraper. Tano Diamond Dusted Membrane Scrapers (DDMS) offer an atraumatic alternative to find the edge of a membrane. The diamond-dusted, soft, silicone tip identifies and lifts the edge of the membrane quickly and easily for your forceps. Available in several styles, including a retractable version, the patented DDMS has become a staple in operating rooms around the world.

DDMS benefits:

- Effectively aids in the atraumatic removal of immature membranes and pigmented cells from the retina
- Assists atraumatic separation of hyaloids from neurosensory retina
- Especially beneficial for relieving traction from the rim of macular holes or retinal breaks in an atraumatic fashion

DDMS™ Diamond Dusted Membrane Scrapers



Tano

- 20.04, 20.04.23 and 20.04.25S
- Available in 20ga, 23ga, and 25ga
- 5/Box



Stiff Tano

- 20.07
- Available in 20ga
- 5/Box



Retractable Tano

- 20.16, 20.16.23 and 20.16.25
- Available in 20ga, 23ga and 25ga
- 20ga — 3/Box,
23ga and 25ga — 5/Box



Adjustable Tano with Pick

- 20.18
- Available in 20ga
- 3/Box



Synergtics® has packaged everything you need for office injections into one complete, sterilized, I-PACK. The I-PACK will save you hours of ordering time, reduce and organize your inventory, and save preparation time since everything you need is conveniently packaged together. The handy plastic tray also serves as a reservoir that you can use for your povidone-iodine solution. Simply pour the solution into the compartment with the speculum, where it can be used for both your speculum and for dipping cotton swabs. The two additional compartments can be used for discarded vials, swabs, and needles, making cleanup quick and easy.

Benefits:

- **SAVE ORDERING TIME:** Affordable, one-stop, quick ordering for all of your office injection requirements.
- **CONTROL INVENTORY:** With only the I-Pack to monitor, keeping track of inventory has never been so easy.
- **SAVE PREPARATION TIME:** No more fumbling around trying to locate misplaced speculums. Everything you need is contained in a ready-to-use, seal-tight, sterilized, package.



Ordering Information

18060 Components

- 1mL Tuberculin Syringes
- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

18061 Components

- 2 1mL Tuberculin Syringes
- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")

18063 Components

- 2 1mL Tuberculin Syringes
- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker
- 1 Disposable Forceps, 4"
- 4 4"x4" Gauze Pads

18064 Components

- 2 1mL Tuberculin Syringes
- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")
- 1 Disposable Forceps, 4"
- 4 4"x4" Gauze Pads

18066 Components

- 2 1mL Tuberculin Syringes
- 1 Pr Diamond Edge High Grip Gloves
- 1 30ga 1/2" Needle
- 1 27ga 1/2" Needle
- 2 18ga 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

18067 Components

- 1 30ga 1/2" Needle
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

18069 Components

- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

18069W Components

- 1 Pr Synthetic Latex-Free Gloves
- 1 30ga 1/2" Needle
- 6 Cotton Swabs 6"
- 1 Wire Eyelid Speculum, Adult
- 1 3.5-4mm Scleral Marker

BIPOLAR [STERILE, 6/BOX]



Bipolar Cautery

- 53.00.23 and 53.00.25
- Available in 23ga and 25ga

KNIVES [STERILE, 5/BOX]



Opremcak BRVO

- 52.00
- Super-thin knife perfect for sheathotomy
- Available in 20ga



Opremcak CRVO

- 52.20
- Ideal for relaxing scleral ring, cribriform plate, and sclera of optic disk during radial optic neurotomy
- Available in 20ga

IRIS RETRACTORS



Reusable

- 40.01-2 and 40.01-6
- 5 retractors/set in case one is lost or damaged
- Autoclavable up to 3 uses
- Available in 2/box or 6/box configurations



Disposable

- 40.02-2 and 40.02-6
- 5 retractors/set in case one is lost or damaged
- Sterile, 2/box or Sterile, 6/box configurations

VISCO DISSECTOR



Hubbard Visco Dissector

- 12.06.25
- Injects viscoelastic for delaminating large membranes from the retina
- Available in 25ga
- Sterile, 6/Box

LENS RINGS



Awh Sutureless Lens Ring

- 41.00
- Used to stabilize a vitrectomy lens
- 3 sizes/pack
- Sterile, 10/Box

DISPOSABLE LENSES [STERILE, 10/BOX]



Flat

- SDLF
- Field of view: 36°
- Image magnification: 1.0x



Bi-Concave

- SDLB
- Field of view: 25°
- Image magnification: 0.8x



Magnifying

- SDLM
- Field of view: 30°
- Image magnification: 1.5x



Widefield

- SDLW
- Field of view: 48°
- Image magnification: 0.5x



30° Prism

- SDL3P
- Field of view: 33° (offset 30°)
- Image magnification: 1.0x



www.synergeticsusa.com