

Exclusively from **Labtician**, the **MST**

# FAB5

**Simplify** the Complicated

Discover the **power of precision**:

Essential tools that help you manage planned and unplanned complications in ophthalmic surgery.



**Packer/Chang IOL Cutters**

Easily cuts any type of foldable IOL.



**Micro-Holding Forceps**

Can firmly grasp the IOL during explantation or repositioning of a dislocated lens.



**Hoffman/Ahmed Horizontal Micro-scissors**

Cuts tissue and sutures within the anterior chamber.



**Ahmed Micro Graspers**

Manipulates tissue and sutures in the anterior chamber.



**Seibel Capsulorhexis Forceps**

Provides superior visibility and control during the capsulorhexis and can save a capsulorhexis that has extended beyond its intended size.



Scan to watch how each of these micro instruments are used in surgery.

For more information or to order, call **1-855-651-4934** or email [orders@labtician.com](mailto:orders@labtician.com)

Exclusively from  
**LABTICIAN**<sup>™</sup>  
Moving Vision Forward

