

F-OCTANE F-DECALIN

NEW! NOW ALSO AVAILABLE
IN READY-TO-USE SYRINGES!



The Surgeon's Third Hand

**Gentle Retinal Unfolding
and Stabilization**

Drainage of Subretinal Fluid

Refloating of luxated lenses

ALSO AVAILABLE!

in 5 and 7 ml vials

COMPARISON OF ULTRAPURE PERFLUOROCARBONS

		F-OCTANE PLUS	F-OCTANE	F-DECALIN PLUS	F-DECALIN
Density [g/cm ³]		1.72	1.76	1.81	1.93
Vapor Pressure [mbar] at 25°C		16.0	18.5	7.0	8.0
Refractive Index at 20°C		1.2770	1.2700	1.3180	1.3110
Surface Tension [mN/m] at 25°C		14.7	14.0	19.4	19.0
Interface Tension [mN/m] at 25°C		54.4	55.0	56.1	57.8
5 ml	Vial	G-80405	G-80305	G-80205	G-80105
	Syringe	–	G-80315	–	G-80115
7 ml	Vial	G-80407	G-80307	G-80207	G-80107
	Syringe	–	G-80317	–	G-80117

ACCESSORIES:



PFCL CANNULA

with dual shaft 24 gauge/0.55 mm
lateral tube: 25 gauge/0.5 mm
G-33055



BACKFLUSH HANDPIECE

with silicone chamber and
luer-lock connector
G-37002



CHANG CANNULA

coaxial, outside diameter of inner
tube: 25 gauge/0.5 mm
20 gauge/0.9 mm
G-33057



SINGLE-USE SILICONE CHAMBER

for backflush handpiece, 10 pcs.,
sterile, **G-34806**



SINGLE-USE BACKFLUSH CANNULA

with silicone brush, 5 pcs., sterile

G-34293 in 20 gauge/0.9 mm
G-34294 in 23 gauge/0.6 mm
G-34297 in 25 gauge/0.5 mm