

LABEDGE™



KNIVES & BLADES PRODUCTS

Our **LABEDGE™** portfolio of knives and blades are offered in a wide variety of sizes to meet all your needs during surgery. **LABEDGE** is designed to offer best-in-class sharpness, precise cutting dimensions, and consistent cutting force through the use of electrochemical sharpening and grindless technology. This ensures the cutting edges of the blades are smooth, absent of residual burrs, and creates a true incision size that can be replicated consistently from blade to blade.

Features Of The LABEDGE™ Single-Use Knives	4
STAB	5
Stab And Paracentesis	5
Controlled Depth	5
CRESCENT	6
Straight	6
Angled	6
PHACO SLIT	7
Single Bevel Angled	7
Double Bevel Angled	7
MVR/SIDEPART	8
Straight	8
Angled	8

Continued next ▶



IMPLANT **8**
 Angled 8

OTHER **9**
 Scleral 9
 Lamellar 9
 Needle 9

MINI BLADES **10**
 Chisel Tip 10
 Spear Tip 10
 Round Tip 10
 Lance Tip 11
 Curved Tip 11
 Round Tip Straight 11
 Round Tip Angled 11


HANDLES **12**
 Cross-cut collet 12
 Round 12
 Hex 13



FEATURES OF THE LABEDGE™ SINGLE-USE KNIVES

The **LABEDGE™** knives have best-in-class sharpness* that reduces the force required for penetration, creating a precise wound architecture and minimizing surface damage.


LABEDGE™ knives are manufactured using electrochemical sharpening and grindless technology to eliminate residual burrs on the cutting edges, producing a true incision size that can be replicated consistently from blade to blade.



For greater incisional accuracy, tip size is indicated in the description and on the handle.

* Testing was performed on a Mark-10 Force Testing gauge and motorized stand that measures force over time.

STAB STAB AND PARACENTESIS



Full Handle Knife


Item No.	Description
LAB961501*	15 °
LAB962201*	22.5 °
LAB963001*	30 °
LAB964501*	45 °

Screw-in Blade

Item No.	Description
LAB961502*	15 °

6 pcs. per package

STAB CONTROLLED DEPTH



Full Handle Knife

Item No.	Description
LAB7514*	15° 3.0 mm depth light blue
LAB7516*	15° 5.0 mm depth Green
LAB7523*	22.5° 4.0 mm depth Grey
LAB7531*	30° 3.5 mm depth Black
LAB7536*	30° 5.0 mm depth Black


Screw-in Blade

Item No.	Description
LAB7511	15° 1.5 mm depth Orange
LAB7513	15° 3.0 mm depth light blue
LAB7515	15° 5.0 mm depth Green
LAB7522*	22.5° 4.0 mm depth Grey
LAB7530*	30° 3.5 mm depth Black

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

CRESCENT STRAIGHT




Full Handle Knife

Item No. Description

LAB960021* 2.0 mm

6 pcs. per package

CRESCENT ANGLED



Full Handle Knife

Item No. Description


LAB940002* 2.5 to 2.0 mm Tapered

LAB960022 2.0 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

PHACO SLIT SINGLE BEVEL ANGLED




Full Handle Knife

Item No.	Description	Item No.	Description
LAB961861*	1.8 mm	LAB962751*	2.75 mm
LAB962061*	2.0 mm	LAB962861*	2.80 mm
LAB962261*	2.2 mm	LAB962851*	2.85 mm
LAB962461*	2.4 mm	LAB962961*	3.0 mm
LAB962561*	2.5 mm	LAB963261*	3.2 mm
LAB962661*	2.65 mm	LAB963561*	3.5 mm

6 pcs. per package

PHACO SLIT DOUBLE BEVEL ANGLED




Full Handle Knife

Item No.	Description	Item No.	Description
LAB961868*	1.8 mm	LAB962758*	2.75 mm
LAB962068*	2.0 mm	LAB962868*	2.80 mm
LAB962268*	2.2 mm	LAB962858*	2.85 mm
LAB962468*	2.4 mm	LAB962968*	3.0 mm
LAB962568*	2.5 mm	LAB963268*	3.2 mm
LAB962668*	2.65 mm		

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


MVR/SIDEPORT STRAIGHT



Item No.	Description	Item No.	Description
LAB965564*	25 ga	LAB965560*	20 ga
LAB965562*	23 ga	LAB965561	19 ga
LAB965563*	21 ga	LAB965550*	1.0 mm incision size

6 pcs. per package


MVR/SIDEPORT ANGLED



Item No.	Description	Item No.	Description
LAB965569*	25 ga	LAB965565*	20 ga
LAB965567*	23 ga	LAB965566*	19 ga
LAB965568*	21 ga	LAB965570*	1.0 mm incision size

6 pcs. per package

IMPLANT ANGLED




Item No.	Description
LAB965261*	5.2 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

OTHER SCLERAL



Full Handle Knife


Item No.	Description
LAB5710*	Standard US style
LAB5710S*	Small European style

Screw-in Blade

Item No.	Description
LAB5700	Standard US style
LAB5700S*	Small European style

6 pcs. per package

OTHER LAMELLAR




Screw-in Blade

Item No.	Description
LAB6600	Angled, bevel up

6 pcs. per package

OTHER NEEDLE




Screw-in Blade

Item No.	Description
LAB5910	3.0 mm
LAB5920*	5.0 mm

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


MINI BLADES CHISEL TIP



Screw-in Blade

Item No.	Description
LAB6100*	Chisel 1.5 mm
12 pcs. per package	


MINI BLADES CHISEL TIP



Screw-in Blade

Item No.	Description
LAB6200*	Chisel 2.7 mm
12 pcs. per package	


MINI BLADES SPEAR TIP



Screw-in Blade

Item No.	Description
LAB6300*	Spear tip
6 pcs. per package	

MINI BLADES ROUND TIP




Screw-in Blade

Item No.	Description
LAB6400	Round tip, sharp one side
12 pcs. per package	

*The lead times to get this product are variable. Please contact your sales representative for more information.

MINI BLADES LANCE TIP




Screw-in Blade

Item No. Description

LAB6500 Lance tip

12 pcs. per package

MINI BLADES CURVED TIP




Screw-in Blade

Item No. Description

LAB6700 Curved tip, sharp one side

12 pcs. per package

MINI BLADES ROUND TIP STRAIGHT




Screw-in Blade

Item No. Description

LAB6900 Round tip, sharp all around, straight

6 pcs. per package

MINI BLADES ROUND TIP ANGLED



Screw-in Blade

Item No. Description

LAB6910* Round tip, sharp all around, angled

6 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

HANDLES CROSS-CUT COLLET



1305

Screw-in Blade

Item No.	Description
LAB1305*	Cross-cut collet
1 pcs. per package	

HANDLES ROUND



1311



1312



1310



1350



1360



1320

Screw-in Blade

Item No.	Description
LAB1311*	5.0 cm round
LAB1312*	7.5 cm round
LAB1310*	10.0 cm round
LAB1350*	13.0 cm round
LAB1360*	13.0 cm round heavy duty
LAB1320*	15.5 cm round
1 pcs. per package	

*The lead times to get this product are variable. Please contact your sales representative for more information.

HANDLES HEX



1330



1340

Screw-in Blade

Item No.	Description
LAB1330*	10.0 cm hex
LAB1340*	15.5 cm hex

1 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.