


LABTICIAN™ LABSPEAR™

Labtician Ophthalmics Lab™SPEAR Catalogue



FEATURES OF THE Lab™SPEAR EYE SPEARS

Lab™SPEAR eye spears are composed of a polyvinyl acetal (PVA) sponge attached to a high-density plastic handle. The sponge is soft and highly absorbent to effectively clean the surgical field while remaining comfortable for the patient. In addition, it can be cut and shaped according to your clinical needs without leaving behind any debris.



HIGH RATE OF WATER ABSORPTION

LAB™SPEAR eye spears can quickly absorb high volumes of liquid, **fully expanding in less than 10 seconds.**

Lab™SPEAR EYE SPEARS

Item No.	Description
LABES1001A 100	With handle, 17 x 10 x 9 mm
100 pcs. per package	

HIGH TENSILE STRENGTH

Research demonstrates that when cutting cellulose sponge eye spears there is a greater chance of micro fragments being left behind.¹ LAB™SPEAR eye spears are manufactured from PVA and can be cut into different shapes without breaking apart or the shedding of fibres thus, **reducing the risk of infection.**¹

SMOOTH AND GENTLE FLUID REMOVAL

The PVA sponge used in LAB™SPEAR eye spears are extremely soft to the touch which allows for **gentle manipulation of delicate tissues** such as the ocular surface.

Labtician Ophthalmics, Inc.
 2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5
 Telephone: 905.829.0055 Toll free: 1.800.265.8391 (Canada only)
 labtician.com

Order Desk & Customer Service
 Toll free: 1.855.651.4934 (Canada & USA only)
 Toll free fax: 1.855.390.9790 (Canada & USA only)
 Email: orders@labtician.com

1. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC1771441/>

© 2024 Labtician Ophthalmics, Inc. All rights reserved.