

 LABTICIAN™

LAB™ TREPHINE

Labtician Ophthalmics Lab™ TREPHINE Catalogue



TREPHINE PRODUCTS

The **Lab™TREPHINE** product line offers a wide selection of corneal trephine blades and donor punches. All trephine products are available in a variety of different sizes and the option to include vacuum functionality, providing you with the ability to secure the trephine firmly in place. In addition, our advanced trephine blade technology is designed specifically to promote reproducible wound architecture while reducing the loss of endothelial cells during transplantation.

LAB™TREPHINE VACUUM TECHNOLOGY	3
RECIPIENT	2
Recipient Advanced Radial Vacuum Trephines	4
DONOR	3
Donor Vacuum Trephine Punches	3
Donor Non-Vacuum Trephine Punches	3
Donor Long Trephines	4
Donor Short Trephine	4
MISCELLANEOUS	5
Dalk Cannula	5
Corneal Block	5



www.labtician.com/docs












LAB™ TREPHINE VACUUM TECHNOLOGY

Stabilization from the vacuum trephines allow the natural contour of the cornea to be maintained during trephination, thereby facilitating a symmetrical and consistent cut each time.

The vacuum trephines and punches are disposable instruments, including a single-hand operated syringe, which create the vacuum seal to hold the donor cornea tissue in place during trephination. This technology simplifies corneal transplants and eliminates the need for additional equipment.

Compact, lightweight design for balance



RECIPIENT ADVANCED RADIAL VACUUM TREPHINES

GUARANTEED PRECISION

The secure vacuum facilitates accurate cuts and reproducible post-operative results.

REDUCED RISK OF SURGICALLY INDUCED ASTIGMATISM

16 radial inking prongs provide symmetrical markings to ensure proper suture placement.

Item No.	Description	Item No.	Description
LAB370600*	6.00 mm	LAB370800	8.00 mm
LAB370650*	6.50 mm	LAB370825*	8.25 mm
LAB370700*	7.00 mm	LAB370850*	8.50 mm
LAB370725*	7.25 mm	LAB370875*	8.75 mm
LAB370750	7.50 mm	LAB370900*	9.00 mm
LAB370775*	7.75 mm		

1 pc. per package



*The lead times to get this product are variable. Please contact your sales representative for more information.

DONOR VACUUM TREPHINE PUNCHES

Vacuum fixation ensures a precise circular cut of the corneoscleral donor button

Item No.	Description	Item No.	Description
LAB360600*	6.00 mm	LAB360800	8.00 mm
LAB360625*	6.25 mm	LAB360825*	8.25 mm
LAB360650*	6.50 mm	LAB360850*	8.50 mm
LAB360675*	6.75 mm	LAB360875*	8.75 mm
LAB360700*	7.00 mm	LAB360900*	9.00 mm
LAB360725*	7.25 mm	LAB360925*	9.25 mm
LAB360750*	7.50 mm	LAB360950*	9.50 mm
LAB360775*	7.75 mm		

1 pcs. per package



DONOR NON-VACUUM TREPHINE PUNCHES

Ergonomic housing supports precise cuts of corneoscleral grafts

Item No.	Description	Item No.	Description
LAB340600*	6.00 mm	LAB340800*	8.00 mm
LAB340625*	6.25 mm	LAB340825*	8.25 mm
LAB340650*	6.50 mm	LAB340850*	8.50 mm
LAB340675*	6.75 mm	LAB340875*	8.75 mm
LAB340700*	7.00 mm	LAB340900*	9.00 mm
LAB340725*	7.25 mm	LAB340925*	9.25 mm
LAB340750*	7.50 mm	LAB340950*	9.50 mm
LAB340775*	7.75 mm		

1 pcs. per package



*The lead times to get this product are variable. Please contact your sales representative for more information.

DONOR LONG TREPHINES

Single-use trephines offer **seamless, precision cutting** and reproducible outcomes
 A wide range of sizes from **6.0 mm to 18.0 mm**
16 mm height

Item No.	Description	Item No.	Description
LAB330250*	2.50 mm	LAB330900*	9.00 mm
LAB330300*	3.00 mm	LAB330950*	9.50 mm
LAB330350*	3.50 mm	LAB331000	10.00 mm
LAB330400*	4.00 mm	LAB331050*	10.50 mm
LAB330500*	5.00 mm	LAB331100*	11.00 mm
LAB330550*	5.50 mm	LAB331150*	11.50 mm
LAB330600	6.00 mm	LAB331200*	12.00 mm
LAB330625*	6.25 mm	LAB331250*	12.50 mm
LAB330650*	6.50 mm	LAB331300*	13.00 mm
LAB330675*	6.75 mm	LAB331350*	13.50 mm
LAB330700	7.00 mm	LAB331400*	14.00 mm
LAB330725	7.25 mm	LAB331450*	14.50 mm
LAB330750	7.50 mm	LAB331500*	15.00 mm
LAB330775	7.75 mm	LAB331550*	15.50 mm
LAB330800	8.00 mm	LAB331600*	16.00 mm
LAB330825	8.25 mm	LAB331650*	16.50 mm
LAB330850	8.50 mm	LAB331700*	17.00 mm
LAB330875	8.75 mm	LAB331800	18.00 mm



1 pcs. per package

DONOR SHORT TREPHINE

Item No.	Description
LAB310800*	8.00 mm




1 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

DALK CANNULA

Item No.	Description
LAB211050*	27G (0.4 mm) Big Bubble, curved, 1 port (Fontana)
10 pcs. per package	



Tapered tip, port on underside

CORNEAL BLOCK

Item No.	Description
LAB330001*	Disposable Corneal Block
1 pc. per package	



*The lead times to get this product are variable. Please contact your sales representative for more information.

Labtician Ophthalmics, Inc.

2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5

Telephone: **905.829.0055** Toll free: **1.800.265.8391** (Canada only)

labtician.com

Order Desk & Customer Service

Toll free: **1.855.651.4934** (Canada & USA only)

Toll free fax: **1.855.390.9790** (Canada & USA only)

Email: **orders@labtician.com**