


LABTM LidLoads

The next generation
of gold eye lid weights
that provide corneal
protection with
enhanced cosmesis

A 3-dimensional solution
for the treatment of
Lagophthalmos

Thickness
Length
Width



PROPORTIONALTM Lid Loads aren't just thinner than standard weights – they were designed and proportioned taking into account all three dimensions (**thickness, length and width**) to fit your patient's needs precisely.

- **Thickness** - thinner in the low weight range for patients with smaller and thinner lids
- **Length** - shorter at the high weight range where managing length is important
- **Width** - rounded perimeter and corners for enhanced cosmesis

Bringing innovation to practice.

PROPORTIONAL[™] Lid Load Gold Eyelid Implants - A unique 3-dimensional design solution

Labtician Standard Gold Weights (mm)			Labtician PROPORTIONAL [™] Gold Weights (mm)		
Product#	Weight	Thickness/Width/Length	Product#	Weight	Thickness/Width/Length
			SGP06	0.6 g	0.55/5.25/12.5
			SGP08	0.8 g	0.575/5.25/15.0
			SGP10	1.0 g	0.60/5.5/16.6
			SGP12	1.2 g	0.675/5.75/16.8
			SGP14	1.4 g	0.70/6.0/17.6
			SGP16	1.6 g	0.725/6.5/17.8
			SGP18	1.8 g	0.75/6.5/19.0
			SGP20	2.0 g	0.78/6.7/20.0
SG22	2.2 g	1.0/5.0/22.0			
SG24	2.4 g	1.0/5.0/23.5			
SG26	2.6 g	1.0/5.0/24.5			
SG28	2.8 g	1.0/5.0/25.5			

Compare with the eye lid weights you are using now

PROPORTIONAL[™]

PROPORTIONAL[™] Lid Loads provide:

1. The corneal protection you want
2. With the enhanced cosmesis your patient expects

Labtician Ophthalmics, Inc.
 2140 Winston Park Drive,
 Unit 6, Oakville, ON,
 Canada L6H 5V5
 Tel: (905) 829-0055
 Fax: (905) 829-0056
 1 800 265 8391
 info@labtician.com
 www.labtician.com



FM72835