


Introducing Malyugin Ring 2.0

Nothing's Changed... Except For Everything



When contemplating what the next version of the Malyugin Ring should provide surgeons and their patients, we established three goals:

1. Provide the same or better safety as surgeons have experienced and have made the Malyugin Ring “classic” the standard of care for pupil management.
2. Make it easier to place and remove from the pupil margin.
3. Allow for entry in 2.0mm incisions and more room in wounds larger than that.

After several years of work, we are proud to introduce the Malyugin Ring 2.0. The second generation of the Malyugin Ring that also can be used in 2.0mm incisions.

Malyugin Ring 2.0 requires no surgical technique changes for you, and we believe that you will find it easier to use, that you will appreciate the extra room it affords, and be glad to know that it has softer compression characteristics (in fact Malyugin Ring 2.0 exerts less than half as much pressure on the iris than the “classic” Malyugin Ring).

Malyugin Ring 2.0

Only from MST

 *Malyugin Ring 2.0*