

# LABTICIAN™ LAB™ LID LOADS

The **Next Generation** of gold eyelid weights that provide corneal protection with enhanced cosmesis

Labtician **Proportional™** Lid Loads aren't just thinner than standard weights — they were designed and proportioned taking into account all three dimensions (**thickness, length and width**) to fit your patient's needs precisely.




**Thickness** - thinner in the low weight range for patients with smaller and thinner lids.



**Length** - shorter at the high weight range where managing length is important



**Width** - rounded perimeter and corners for enhanced cosmesis



Unique 3-dimensional design solution



A unique 3-dimensional design solution



### Labtician Proportional™ Gold Weights (mm)

Product	Weight	Thickness	Width	Length
SGP06	0.6 g	0.55 mm	5.25 mm	12.5 mm
SGP08	0.8 g	0.575 mm	5.25 mm	15.0 mm
SGP10	1.0 g	0.60 mm	5.5 mm	16.6 mm
SGP12	1.2 g	0.675 mm	5.75 mm	16.8 mm
SGP14	1.4 g	0.70 mm	6.0 mm	17.6 mm
SGP16	1.6 g	0.725 mm	6.5 mm	17.8 mm
SGP18	1.8 g	0.75 mm	6.5 mm	19.0 mm
SGP20	2.0 g	0.78 mm	6.7 mm	20.0 mm

Compare with the eyelid weights you are using now

Proportional™ Lid Loads provide:

1. The corneal protection you want
2. With the enhanced cosmesis your patient expects



intertek

Labtician Ophthalmics, Inc.  
 2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5  
 Telephone: 905.829.0055 Toll free: 1.800.265.8391 (Canada only)  
 labtician.com

Order Desk & Customer Service  
 Toll free: 1.855.651.4934 (Canada & USA only)  
 Toll free fax: 1.855.390.9790 (Canada & USA only)  
 Email: orders@labtician.com