

 LABTICIAN™

LAB™CANNULA

Labtician Ophthalmics Lab™CANNULA Catalogue



CANNULA PRODUCTS

We offer a diverse portfolio of cannulas of the highest quality, ranging from 19g to 30g to meet your various needs. **Lab™ CANNULA** are manufactured in North America, where a high level of consistency and product quality is maintained to ensure patient safety and promote successful outcomes.

FEATURES OF THE LAB™ CANNULA SINGLE-USE CANNULAS	1
ANESTHESIA	2
Retrobulbar, Atkinson Tip	2
Retrobulbar, Sharp Tip	2
Peribulbar, Atkinson Tip	2
Sub-Tenon	2
ANTERIOR CHAMBER IRRIGATION	3
Angled Tip, 4mm (Rycroft)	3
Angled Tip, 6mm (Rycroft)	3
Angled Tip, 7mm (Rycroft)	3
Angled Tip, 8mm (Lasik)	3
Angled Tip, 10mm (Lasik)	4
CAPSULE POLISHING	5
Silicone Tip, Angled 7mm	5
Posterior, Angled 9mm, Top Port	5
Posterior, Angled 7mm, Open Tip	5
Olive Tip	5

Continued next ▶

CAPSULOTOMY	6
Formed	6
Angled 45°	6
Reverse Bend	6
Short	6
Side-Cut	7
Push-Pull, Formed	7
HYDRODISSECTION	8
Flat Tip, Angled 7mm	8
Flat Tip, Angled 8mm	8
Horizontal Flat Tip, Curved	8
Vertical Flat Tip, Curved	8
J-Shaped, Vertical Flat Tip	9
Angled Standard Tip (Chang)	9
Cortical Cleaving, Flat Tip, Angled 11mm	9
Cortical Cleaving, Round Tip, Angled 11mm	9
Cortical Cleaving, Round Tip, Angled 12mm	10
Tapered, Angled 1mm	10
IRRIGATION AND AIR INJECTION	11
Viscoelastic, Angled	11
Viscoelastic, Formed	11
Washout, Angled 8mm (Bishop-Harmon)	11
Olive Tip, Curved	11
NUCLEUS REMOVAL	12
Vectis, 1 Port	12
IRRIGATION AND ASPIRATION	12
Simcoe, Standard	12
ANTERIOR CHAMBER MAINTENANCE	12
Ribbed Open Tip	12



www.labtician.com/docs




FEATURES OF THE LAB[™]CANNULA SINGLE-USE CANNULAS

CANNULAS THAT MEET YOUR DIVERSE SURGICAL NEEDS

- Labtician provides a wide range of **innovative cannula designs** with the highest quality.
- Our cannulas are manufactured in North America where a high level of **consistency and product quality** is maintained.
- Single-Use surgical instruments are highly **recommended by the ASCRS** (American Society for Cataract & Refractive Surgeons)¹. Instruments such as cannulas contain narrow tubing and are immensely difficult to clean. Any errors in the cleaning process can increase the risk of potential infections including TASS (toxic anterior segment syndrome). By utilizing single-use cannulas you are reducing this risk and increasing the probability of a successful patient outcome.


Labtician Consumables hubs are colour-coded in accordance with ISO 6009 standards, with the exception of 30G (clear).



Colour	Gauge	Description	Colour	Gauge	Description
CLEAR	30G	0.30 mm	BLUE	23G	0.60 mm
GREY	27G	0.40 mm	BLACK	22G	0.70 mm
BROWN	26G	0.45 mm	GREEN	21G	0.80 mm
ORANGE	25G	0.50 mm	YELLOW	20G	0.90 mm
PURPLE	24G	0.55 mm	CREAM	19G	1.10 mm

¹ http://tassregistry.org/TASS_guidelines-CBC.pdf


ANESTHESIA RETROBULBAR, ATKINSON TIP



Item No.	Description	
LAB111638*	23G x 1 1/2 in.	0.6 x 38 mm
LAB111637	25G x 1 1/2 in.	0.5 x 38 mm

10 pcs. per package


ANESTHESIA RETROBULBAR, SHARP TIP



Item No.	Description	
LAB111278*	25G x 1 in.	0.5 x 25 mm
LAB111282*	30G x 1/2 in.	0.3 x 12.5 mm

10 pcs. per package


ANESTHESIA PERIBULBAR, ATKINSON TIP



Item No.	Description	
LAB115144*	23G x 7/8 in.	0.6 x 22 mm
LAB115027*	25G x 7/8 in.	0.5 x 22 mm
LAB111642*	25G x 1 1/4 in.	0.5 x 32 mm

10 pcs. per package

ANESTHESIA SUB-TENON




Item No.	Description	
LAB125176	19G x 7/8 in.	1.1 x 22mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


ANTERIOR CHAMBER IRRIGATION ANGLED TIP, 4 MM (RYCROFT)



Item No.	Description	
LAB121281*	20G	0.9 mm
LAB121639*	25G	0.5 mm
LAB101280	27G	0.4 mm
LAB121273	30G	0.3 mm

10 pcs. per package


ANTERIOR CHAMBER IRRIGATION ANGLED TIP, 6 MM (RYCROFT)



Item No.	Description	
LAB121287*	20G	0.9 mm
LAB121286*	27G	0.4 mm
LAB121276*	30G	0.3 mm

10 pcs. per package


ANTERIOR CHAMBER IRRIGATION ANGLED TIP, 7 MM (RYCROFT)



Item No.	Description	
LAB121637*	25G	0.5 mm
LAB121290*	27G	0.4 mm
LAB121277*	30G	0.3 mm

10 pcs. per package

ANTERIOR CHAMBER IRRIGATION ANGLED TIP, 8MM (LASIK)




Item No.	Description	
LAB121293*	27G	0.4 mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.

ANTERIOR CHAMBER IRRIGATION ANGLED TIP, 10MM (LASIK)




Item No.	Description	
LAB121280*	25G	0.5 mm
LAB121274*	27G	0.4 mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


CAPSULE POLISHING SILICONE TIP, ANGLED 7MM



Item No.	Description	
LAB175041*	23G	0.6 mm
LAB175005*	27G	0.4 mm

10 pcs. per package


CAPSULE POLISHING ANGLED 9MM, TOP PORT



Item No.	Description	
LAB171604*	21G	0.8 mm

10 pcs. per package


CAPSULE POLISHING ANGLED 7MM, OPEN TIP



Item No.	Description	
LAB171614*	21G	0.8 mm

10 pcs. per package

CAPSULE POLISHING OLIVE TIP




Item No.	Description	
LAB171589*	23G	0.6 mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


CAPSULOTOMY FORMED



Item No.	Description	
LAB131599*	23G x 5/8 in.	0.6 x 16 mm
LAB131610	25G x 5/8 in.	0.5 x 16 mm
LAB131618	27G x 5/8 in.	0.4 x 16 mm

10 pcs. per package


CAPSULOTOMY ANGLED 45°



Item No.	Description	
LAB135049*	25G x 5/8 in.	0.5 x 16 mm
LAB135093	27G x 5/8 in.	0.4 x 16 mm

10 pcs. per package


CAPSULOTOMY REVERSE BEND



Item No.	Description	
LAB135115	25G x 5/8 in.	0.5 x 16 mm

10 pcs. per package

CAPSULOTOMY SHORT




Item No.	Description	
LAB135004*	25G x 5/8 in.	0.5 x 16 mm
LAB135008*	27G x 5/8 in.	0.4 x 16 mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


CAPSULOTOMY SIDE-CUT



Item No.	Description	
LAB131602*	25G x 5/8 in.	0.5 x 16 mm
LAB131609*	27G x 5/8 in.	0.4 x 16 mm

10 pcs. per package

CAPSULOTOMY PUSH-PULL, FORMED




Item No.	Description	
LAB135042*	27G x 5/8 in.	0.4 x 16 mm

10 pcs. per package

*The lead times to get this product are variable. Please contact your sales representative for more information.


HYDRODISSECTION FLAT TIP, ANGLED 7MM



Item No.	Description	
LAB125024*	23G	0.6 mm
LAB125036*	25G	0.5 mm
LAB125162*	27G	0.4 mm

10 pcs. per package


HYDRODISSECTION FLAT TIP, ANGLED 8MM



Item No.	Description	
LAB125037	25G	0.5 mm
LAB125158*	27G	0.4 mm
LAB121592*	30G	0.3 mm

10 pcs. per package


HYDRODISSECTION HORIZONTAL FLAT TIP, CURVED



Item No.	Description	
LAB125155*	27G	0.4 mm

10 pcs. per package

HYDRODISSECTION VERTICAL FLAT TIP, CURVED




Item No.	Description	
LAB125099*	27G	0.4 mm

10 pcs. per package


*The lead times to get this product are variable. Please contact your sales representative for more information.

HYDRODISSECTION J-SHAPED, VERTICAL FLAT TIP




Item No.	Description		
LAB125029*	27G	0.4 mm	
10 pcs. per package			

HYDRODISSECTION ANGLED STANDARD TIP (CHANG)




Item No.	Description		
LAB121198	27G	0.4 mm	1.5mm tip
10 pcs. per package			

HYDRODISSECTION CORTICAL CLEAVING, FLAT TIP, ANGLED 11MM



Item No.	Description		
LAB105158*	27G	0.4 mm	
10 pcs. per package			


HYDRODISSECTION CORTICAL CLEAVING, ROUND TIP, ANGLED 11MM



Item No.	Description		
LAB125157*	27G	0.4 mm	
10 pcs. per package			


*The lead times to get this product are variable. Please contact your sales representative for more information.

HYDRODISSECTION CORTICAL CLEAVING, ROUND TIP, ANGLED 12MM



Item No.	Description	
LAB125197	30G	0.3 mm
10 pcs. per package		


HYDRODISSECTION TAPERED, ANGLED 1MM



Item No.	Description	
LAB125159	25G - 30G	0.5 - 0.3 mm
10 pcs. per package		

*The lead times to get this product are variable. Please contact your sales representative for more information.


IRRIGATION AND AIR INJECTION VISCOELASTIC, ANGLED



Item No.	Description	
LAB125009*	23G x 7/8 in.	0.6 x 22 mm
LAB125026	27G x 7/8 in.	0.4 x 22 mm

10 pcs. per package


IRRIGATION AND AIR INJECTION VISCOELASTIC, FORMED



Item No.	Description	
LAB125007	27G x 7/8 in.	0.4 x 22 mm

10 pcs. per package


IRRIGATION AND AIR INJECTION WASHOUT, ANGLED 8MM (BISHOP-HARMON)



Item No.	Description	
LAB165040	19G x 5/16 in.	1.1 x 8 mm

10 pcs. per package

IRRIGATION AND AIR INJECTION OLIVE TIP, CURVED




Item No.	Description	
LAB165147*	23G x 7/8 in.	0.6 x 22 mm

10 pcs. per package


*The lead times to get this product are variable. Please contact your sales representative for more information.

NUCLEUS REMOVAL VECTIS, 1 PORT




Item No.	Description
LAB141619*	24G 0.55 mm
10 pcs. per package	

IRRIGATION AND ASPIRATION SIMCOE, STANDARD



Item No.	Description
LAB161571	23G 0.6 mm disposable
LAB161573*	23G 0.6 mm reusable
10 pcs. per package - disposable	
10 pcs. per package - reusable	

ANTERIOR CHAMBER MAINTENANCE RIBBED OPEN TIP



Item No.	Description
LAB125061*	20G 0.9 mm disposable
LAB125062*	23G 0.6 mm disposable
3 pcs. per package	

*The lead times to get this product are variable. Please contact your sales representative for more information.

Labtician Ophthalmics, Inc.

2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5

Telephone: **905.829.0055** Toll free: **1.800.265.8391** (Canada only)

labtician.com

Order Desk & Customer Service

Toll free: **1.855.651.4934** (Canada & USA only)

Toll free fax: **1.855.390.9790** (Canada & USA only)

Email: **orders@labtician.com**